


mis® | MGUIDE

FOR C1 & V3 IMPLANTS




STEP-BY-STEP GUIDED SURGICAL PROCEDURE



LOGICAL LAYOUT


1 Example: C1 Implant Ø4.20 / 11.5L; **Soft Bone**



Ø4.20 / 11.5L




② Example: Implant C1 Ø4.20 / 11.5L ; **Hard Bone**



Ø4.20 / 11.5L



PRE-SURGICAL STEPS



VERIFICATION

The package includes:

- A surgical template
- Documentation; including information specific to each planned implant.

KITS REQUIRED

The kit selection is marked on the MGUIDE box label.

- Prior to surgery: Ensure that surgical template, plan and documentation are all made according to the doctor's specifications, and for the relevant patient.




DISINFECTION

The MGUIDE template is shipped non-sterile. Therefore, the template must undergo disinfection prior to use.

DISINFECTION: Completely immerse in a 0.2% Chlorhexidine solution for 10 minutes at room temperature prior to surgery.

- **WARNING!** Do not autoclave. Steam sterilization will deform the template.



INITIAL TRY-IN

It is essential to try-in the template in the patient's mouth, prior to surgery. Correct seating and stability of the template must be confirmed, as well as sufficient space for surgical tools.

In order to avoid incorrect seating due to patient's anatomy change, MGUIDE Template should be used within 3 months from CBCT scan date. No changes shall be made to the oral cavity unless discussed with the MCENTER.

- In rare cases, minor adjustments may be required.



PRECAUTIONS

GENERAL

- All MGUIDE drills and instruments are for use ONLY with the MGUIDE surgical template.
- Metal sleeves must be firmly attached to the template.
- Inspect all instruments prior to each surgery and replace if broken or dull.
- Ensure cooling of cutting instruments with sterile saline solution.
- Tissue punch is NOT equipped with built-in stoppers.
- The MGUIDE Kit to be used, is marked on the MGUIDE Box label.

HANDLING

- Hold the template firmly while drilling.
- Drills and tools MUST engage the sleeve before contra-angle is activated.
- Avoid lateral pressure on the instruments, as it may result in a shift in template position, detachment of sleeves from the template or damage to instruments.
- Use an 'in-out' motion while drilling, slowly inserting the drill until the built-in stopper touches the sleeve.
- Do not over-tighten implant insertion tools and fixation pins. This may result in a shift in template position or damage to the template.



Drills and tools **MUST** engage the sleeve before contra-angle is activated.

FIXATION PINS

(When applicable)



MG-DFP20

MGUIDE drill for fixation pin, Ø2mm



MG-FP020

MGUIDE fixation pin, Ø2mm

FIXATION PINS

MGUIDE fixation pins are recommended for use in fully edentulous cases or if template stability cannot be guaranteed.



- Fixation pins may **ONLY** be used when included in the MGUIDE surgical plan and when pin location is guided by the template.
- Template **MUST** be verified in position, and held firmly prior to drilling.
- Use the MG-DFP20 drill **ONLY!** Drill until stopper touches the sleeve.

TISSUE PUNCH



TISSUE PUNCH

The tissue punch creates a round cut beneath the sleeve. This marks the implant position.




TISSUE REMOVAL

Remove the template, then manually remove punched gingiva.

- Leave at least 2mm of attached gingiva around each implant site.
- Tissue punch tools DO NOT have built-in stoppers.


BONE MILL

(When applicable)



BONE MILL


The bone mill is designed to flatten the alveolar ridge, when needed, prior to drilling.



BONE MILL USE

A flat surface allows for a better approach for the starter drill, therefore increasing the accuracy for the rest of the drilling sequence.

- Use of a bone mill should be part of the planning stage.




BONE MILL STOPPER

A built-in stopper is used for depth control. Bone mill drills have a built-in stopper.

BONE ANCHOR OSTEOTOMY




(When applicable)



STARTER DRILL

Anchor screws are used to vertically secure the template into an osteotomy created by the starter drill MG-D0624.



ANCHORING

Anchor screws should be placed manually. If needed, secure screws using a ratchet until stopper touches the sleeve.

Ø5.5mm sleeve	Ø4mm sleeve
CG-TAS55	CG-NTAS0

- Do not over-tighten screws, as this may cause damage to the template.

PROCEDURE



Drilling Speed (RPM)	400-600	400-600	Soft Bone (type 3&4)	Hard Bone (type 1&2)	15-25
Diameter	Ø2.40	Ø2.40	Ø2.30 Ø3.25	Ø2.30 Ø3.25	Ø3.30

Drilling Speed (RPM)	400-600	400-600	400-600	200-400	15-25
Diameter	Ø2.40	Ø2.40	Ø2.30 Ø3.25	Ø2.90 Ø3.85	Ø3.75

- Do not use the last drill for bone types 3&4.
- The drilling sequence is demonstrated by a 11.5mm C1 implant.
- Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.

PROCEDURE



Drilling Speed (RPM)		Soft Bone (type 3&4)						Hard Bone (type 1&2)	
Diameter		Ø2.40	Ø2.40	Ø2.30 Ø3.25	Ø2.90 Ø3.85	Ø3.40 Ø4.25	Ø4.20	Ø2.40	Ø2.40
Ø4.20/ Ø4.30	MGUIDE Template								

Drilling Speed (RPM)		400-600		400-600		400-600		400-600	
Diameter		Ø2.40	Ø2.40	Ø2.30 Ø3.25	Ø2.90 Ø3.85	Ø3.40 Ø4.25	Ø3.90 Ø4.95	Ø5	
Ø5	MGUIDE Template								

- Do not use the last drill for bone types 3&4.
- The drilling sequence is demonstrated by a 11.5mm C1 implant.
- Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.

DIRECT RATCHET INSERTION TOOL




(When applicable)



IMPLANT INSERTION

The direct insertion tool should be attached to the implant manually prior to implant placement. It should remain connected to the implant until the implant placement procedure has been completed, and then removed manually.

- This option is valid only if the implant has significant primary stability.



STABILIZATION

The direct insertion tool may be used to stabilize the MGUIDE.

- Do not over-tighten the direct insertion tool, as this may cause damage to the template.

IMPLANT INSERTION OPTIONS



BY MOTOR

Recommended for initial implant placement.

Ø5.5mm sleeve	Ø4mm sleeve
CG-GMS10	CG-NMS10
CG-GMN10	CG-NMN10
CG-GMW10	VG-NMN10
VG-GMN10	



BY RATCHET

Orientation adjustment may be achieved only after template removal.

Ø5.5mm sleeve	Ø4mm sleeve
CG-GRS10	CG-NRS10
CG-GRN10	CG-NRN10
CG-GRW10	VG-NRS10
VG-GRS10	VG-NRN10
VG-GRN10	



SAFETY TOOL



BY RATCHET - DIRECT

In cases where additional stability is required for the template.


Ø5.5mm sleeve	Ø4mm sleeve
CG-GRS01	CG-NRS01
CG-GRN01	CG-NRN01
CG-GRW01	VG-NRN01
VG-GRN01	



The guided drill length gauge (MG-DLG55), verifies drill length and may be used before, during and after surgery.

- Taking measurements: Place drill stopper in contact with the gauge. Measure to the drill tip.

DIRECT RATCHET DISASSEMBLY



After each use, the ratchet wrench's adapter should be removed and direct insertion tools should be disassembled prior to cleaning. Reassembly prior to sterilization is required.

EXTRACTION PROCEDURE

A	Assembled
B	Loosen the screw (counter clockwise)
C	Removal of assembly

- When assembling, tighten screw by hand only.



- MGUIDE Template should be used within 3 months from CBCT scan date.
- For cleaning and sterilization instructions, please refer to the 'Cleaning and Maintenance Instructions for Surgical Instruments' included.

Key to codes used:

	Batch code
	Catalog number
	Manufacturer
	Non-sterile



MIS Implants Technologies Ltd.
P.O.Box 7, Bar Lev Industrial Park,
2015600, Israel
Website: www.mis-implants.com



MIS Implants Technologies GmbH
Simeonskarre 2, 32423 Minden, Germany
Tel: +49 571-972-7620
Email: service@mis-implants.de

MGUIDE KITS FOR CONICAL CONNECTION IMPLANTS

MGUIDE Tools Kit for Conical Connection Implants



MGUIDE Drills Kit for Conical Connection Implants





MIS | **M**GUIDE

The MIS Quality System complies with international quality standards: ISO 13485:2016 – Quality Management System for Medical Devices, ISO 9001:2008 – Quality Management System and CE Medical Device Directive 93/42/EEC.