

mis[®]

C1 CASES 2013

CASE 1

P.
CARDELLI
ITALY

C1 AND BONDBONE IN THE ESHTETIC ZONE

Paolo Cardelli

C1 AND BONDBONE

INITIAL



C1 AND BONDBONE

INITIAL



SEVERE DEEP BITE

C1 AND BONDBONE

INITIAL



C1 AND BOND BONE

INITIAL



C1 AND BONDBONE

CBCT TEMPLATE




C1 AND BOND BONE

CBCT



BUCCAL BONE PLATE: $< 1 \text{ mm}$



CEMENTED
RESTORATION

C1 AND BOND BONE

EVALUATION

"Type 2" implant placement

Advantages:

Soft and hard tissue stability

Preservation of a buccal bony contour

Disadvantages:

2-3 surgical procedures: extraction-implant-exposure?

provisional management (*occlusion?*)

Chen ST, Buser D.

Clinical and esthetic outcomes of implants placed in postextraction sites.

Int J Oral Maxillofac Implants. 2009;24 Suppl:186-217.

Buser D, Chappuis V, Bornstein MM, Wittneben JG, Frei M, Belser UC.

Long-Term Stability of Contour Augmentation With Early Implant Placement Following Single Tooth Extraction in the Esthetic Zone. A Prospective, Cross-Sectional Study in 41 Patients With a 5- to 9-Year Follow-Up.

J Periodontol. 2013 Jan 24. [Epub ahead of print]

C1 AND BOND BONE

EVALUATION

Type 1 vs Type 2 implant placement

*Both treatment approaches appear to be appropriate,
with the preferred treatment based on factors
other than resultant soft tissue changes.*

van Kesteren CJ, Schoolfield J, West J, Oates T.

A prospective randomized clinical study of changes in soft tissue position following immediate and delayed implant placement.

Int J Oral Maxillofac Implants. 2010 May-Jun;25(3):562-70.

C1 AND BOND BONE

EVALUATION

Type 1

Predictable technique

Treatment of choice in cases of single anterior tooth

Correct positioning of the implants

Maintaining the original condition of both bone and soft tissues around the tooth

Malchiodi L, Cucchi A, Ghensi P, Nocini PF. Evaluation of the Esthetic Results of 64 Nonfunctional Immediately Loaded Postextraction Implants in the Maxilla: Correlation between Interproximal Alveolar Crest and Soft Tissues at 3 Years of Follow-Up. Clin Implant Dent Relat Res. 2013 Feb;15(1):130-142.

C1 AND BOND BONE

EVALUATION

Ridge alteration following tooth extraction:

*Phase 1: Bundle bone --> woven bone substitution
(mainly buccal wall) --> vertical reduction*

*Phase 2: resorption from the outer surfaces of
both bone walls.*



Esthetic outcome?

Araújo MG, Lindhe J.

Dimensional ridge alterations following tooth extraction. An experimental study in the dog.

J Clin Periodontol. 2005 Feb;32(2):212-8.

C1 AND BOND BONE

PLANNING

Buccal plate augmentation

Even subtle postextraction buccal plate resorption may have significant clinical effects, particularly in the esthetic zone. Buccal plate augmentation consists of placement of bone graft material over an intact buccal plate, underneath the soft tissues in a surgically created pouch with an aim to maintain or augment the soft tissue esthetics of the region.

Caiazzo A, Brugnami F, Mehra P.

Buccal plate augmentation: a new alternative to socket preservation.

J Oral Maxillofac Surg. 2010 Oct;68(10):2503-6. doi: 10.1016/j.joms.

2010.05.044.

C1 AND BOND BONE

PLANNING

Single flap approach

"Type I" implant placement

Buccal plate augmentation (composite)

Immediate provisional

Trombelli L, Farina R, Franceschetti G, Calura G.

Single-flap approach with buccal access in periodontal reconstructive procedures.

J Periodontol. 2009 Feb;80(2):353-60. doi: 10.1902/jop.2009.080420 .

Caiazzo A, Brugnami F, Mehra P.

Buccal plate augmentation: a new alternative to socket preservation.

J Oral Maxillofac Surg. 2010 Oct;68(10):2503-6. doi: 10.1016/j.joms.2010.05.044.

Casap N, Zeltser C, Wexler A, Tarazi E, Zeltser R.

Immediate placement of dental implants into debrided infected dentoalveolar sockets.

J Oral Maxillofac Surg. 2007 Mar;65(3):384-92.

C1 AND BOND BONE SURGERY

SINGLE FLAP APPROACH



ATRAUMATIC EXTRACTION



C1 AND BOND BONE SURGERY


2.4 TWIST DRILL




OT4 ON THE PALATAL WALL



OT4 WITHOUT TEMPLATE




3.0 TWIST DRILL



C1 AND BOND BONE SURGERY


COUNTERSINK




C1 3.75x10



FINAL POSITION




1 MM SUBCRESTAL



C1 AND BOND BONE



IMMEDIATE PROVISIONAL

IMMEDIATE PROVISIONAL FROM SURGICAL TEMPLATE





C1 AND BONDBONE SURGERY

COMPOSITE GRAFT: BONDBONE (50%) AND β -TCP/HA (50%)



C1 AND BONDBONE SURGERY

COMPOSITE GRAFT: BONDBONE (50%) AND β -TCP/HA (50%)
WITHOUT MEMBRANE



C1 AND BONDBONE SURGERY

PROVISIONAL PLACEMENT + BUCCA GAP FILLING (β -TCP/HA)



SUTURE PGA 6-0



PALATAL ASPECT



C1 AND BONDBONE

HEALING

2 WEEKS



1 MONTH



2 MONTHS



2.5 MONTHS



C1 AND BONDBONE

IMPRESSION (3 MONTHS)




ADEQUATE BUCCAL CONTOUR




C1 AND BONDBONE

IMPRESSION (3 MONTHS)



EMERGENCE PROFILE ACQUISITION



C1 AND BONDBONE

OCCUSAL
RELATIONSHIP



C1 BASE



“VENEER PREPARED”
LITHIUM DISILICATE
ABUTMENT



C1 AND BONDBONE

ABUTMENT TRY-IN




C1 AND BONDBONE

ABUTMENT FINAL SEATING




C1 AND BONDBONE

RESTORATION DELIVERY



C1 AND BONDBONE

FINAL OUTCOME




C1 AND BONDBONE

FINAL OUTCOME




C1 AND BONDBONE

SOFT TISSUE CONTOUR COMPARISON, RIGHT-LEFT



C1 AND BONDBONE

SOFT TISSUE CONTOUR COMPARISON, PRE AND POST-OP



CASE 2

J.ALVAREZ CANTONI ARGENTINA

FUNCTION & ESTHETIC

AN INTERDISCIPLINARY APPROACH

Topic: "The combined use of MIS' dental implants and BONDBONE for immediate procedures in the esthetic zone".

Author: Joile Alvarez Cantoni

ARGENTINA



MIS CLINICAL CASES COMPETITION



Topic: "The combined use of MIS' dental implants and BONDBONE for immediate procedures in the esthetic zone".

Author: Joile Alvarez Cantoni

ARGENTINA

PAGE 1

FUNCTION & ESTHETIC

AN INTERDISCIPLINARY APPROACH

INTRODUCTION

The placement of dental implants in the esthetic zone is a real challenge for doctors because of patients demanding esthetics and difficult anatomic pre-existing terrain.

Potential causes of esthetics failures, pre-operative analysis, ideal implant 3D position and restorative aspects will be discussed in the beginning of this clinical case presentation.

The esthetic zone can be defined as any area to be fixed that is visible in the patient's full smile. An esthetic implant prosthesis is one that resembles a tooth in all aspects. The exact location in which the implant is placed is of extreme importance, it should have the correct location in all three dimensions: apicocoronally, mesiodistally and faciolingually. Any diversion from these positions will have a negative effect in the final restoration.



The esthetic zone

Over the past 15 years, dental esthetics has been an important issue in implant dentistry. In the esthetic zone, unsatisfactory treatment results can lead to devastating clinical situations that can only be reestablished with the removal of the implant and the future surgical ridge and soft tissue augmentation.

In the kick-off of esthetic area implant therapy we should start understanding the patient's desires. In most cases the patient demands an esthetic tooth replacement providing a beautiful smile. It is our responsibility to have the knowledge of all treatment possibilities. Nowadays, implant supported restorations mostly represent the best solution.

Hard tissue deficiencies mostly often need guided bone regeneration to allow the three-dimensional correct implant placement, even if we decide to use a simultaneous or staged approach.

To successfully meet the outcomes of esthetic implant therapy an interdisciplinary team approach is critical and highly recommended.


PAGE 2

3-DIMENSIONAL POSITION OF THE IMPLANT


Placing the implants in a proper 3-dimensional position is a key to an esthetic outcome regardless of the implant system used.

Two anatomical structures have great importance: the bone height of the alveolar crest in the interproximal areas and the height and thickness of the facial bone wall, where the interproximal crest is responsible for the absence or presence of peri-implant papillae.


Several surgical techniques have been presented in the past 15 years to improve bone defects at the facial aspects of implant sites, such as onlay grafting, GRB using barrier membranes, a combination of block bone grafts and barrier membranes, and distraction osteogenesis.



Correct implant position in the orofacial dimension, the implant shoulder is positioned 2 mm palatal to the buccal plate.



Correct implant position in the mesiodistal dimension, the implant shoulder should be positioned 1.5 mm off the neighboring tooth.




Correct implant position in the apicocoronal dimension, the implant shoulder should be positioned 1.5-2 mm apical to the Cement Enamel Junction (CEJ) of the contralateral tooth.

PAGE 3

3-DIMENSIONAL POSITION OF THE IMPLANT CAUSES OF ESTHETIC IMPLANT FAILURE

Iatrogenic factors: failures can be caused by inappropriate implant positioning or implant selection.



Compromised esthetic result in an adult patient. Clinical situation 1 year following implant restoration. The periapical radiograph clearly shows the cause of the esthetic failure: the implant shoulder was positioned too far apically and near the neighboring tooth, which led to the resorption of the buccal plate.



Disastrous result in the esthetic zone, clinical status 6 month following implant placement. The implant shoulders were positioned too facially and close to each other.



Compromised esthetic result in a young female patient. Clinical situation 2 years following implant restoration. The incorrect selection of a wide diameter implant evidences the cause of the esthetic failure. The periapical radiograph clearly shows the resorption of the buccal plate and the interproximal bone.

PAGE 4

CLINICAL TREATMENT PLANNING

CASE PRESENTATION

Age at initial presentation: 20 years

Initial presentation: October 2011

Active treatment completed: February 2013

INTRODUCTION & BACKGROUND

The patient is a 20-year-old laws student and was initially seen by his general dentist with the goal of a simple correction of her composite veneers. She was encouraged to consider an interdisciplinary approach. It was then when she became aware of the need to improve her dental health and her dental esthetics for professional and social reasons, she was willing to accept a plan that would address all of her biological, functional and esthetics needs.

MEDICAL HISTORY

The patient was in excellent health

DIAGNOSTIC FINDINGS

Esthetics analysis: there was no significant alterations in the facial profile.

Incisal plane: convex.

Incisal profile: right maxillary incisor was markedly retrusive, both central incisors were in a protrusive position and the left lateral incisor was in a correct situation.

Incisal length: right lateral incisor 10 mm, right central incisor 10 mm, left central incisor 11.5 mm and left lateral incisor 9.5 mm.

Tooth proportion: non harmonious proportions between the maxillary incisors

Gingival plane: altered gingival plane, showing an asymmetrical appearance in the six anterior maxillary dental elements.

Intraoral dental findings:

- Tooth # 1.7-1.6-2.7-2.6-3.7-3.6-4.7-4.6 Relapsed composite fillings.

- Tooth # 1.2 all-metal ceramic restoration.

- Tooth # 1.1-2.1-2.2 composite veneer.

PAGE 5

CLINICAL TREATMENT PLANNING

PRETREATMENT



Auspicious facial analysis - Favorable low smile line - Great lip support



Incorrect gingival plane - Lack of keratinized gingiva - Inadequate restorations



Incorrect gingival plane - 5 mm probing between #2.1-2.2

PAGE 6

CLINICAL TREATMENT PLANNING

Intraoral periodontal findings:

- Very good plaque control.
- Thin scalloped marginal periodontium
- Chronic gingival inflammation around the right lateral incisor due to marginal overhang
- Lack of keratinized gingiva on the facial aspect of teeth # 1.2-1.1-2.1-2.2.
- Probing depths of teeth # 1.2-1.1 within 3 mm
- Probing depth of tooth # 2.1: Facial 4.5 mm. Distal 6 mm.
- Probing depth of tooth # 2.2: Facial 3 mm. Mesial 5 mm.

Radiographic findings:

- Tooth # 1.2 endodontic failure, lack of gutta-percha condensation. Cast post and core.
- Tooth # 1.1 endodontic failure showing apical radiolucency.
- Tooth # 2.1 endodontic failure showing external resorption.
- Tooth 2.2 endodontic failure showing apical radiolucency.




The periapical radiographs show the endodontic failure of the four maxillary anterior teeth, apical radiolucency and external resorption of number 2.1. It can also be seen the lack of interproximal bone between # 2.1 and 2.2 delivering a critical resolution in the GBR augmentation and implant placement.

PAGE 7

CLINICAL TREATMENT PLANNING

3D ct scans findings:

- Tooth # 1.2-2.1 apical radiolucency.
- Tooth # 2.1 absence of bone buccal plate, extreme protrusive position. Less than 1 mm of interproximal bone in its distal aspect.
- Tooth # 2.1 apical radiolucency and external resorption.
- Poor amount of apical bone on both central incisors.



Endodontic failure - Apical radiolucency - External resorption

PAGE 8

CLINICAL TREATMENT PLANNING

PROPOSED TREATMENT PLANS

REVIEW OF TREATMENT GOALS

1. Harmonize all tooth proportions.
2. Enhance the biological and functional missituation of # 2.1.
3. Return the correct integration of the esthetic zone with the patient's smile.
4. Eliminate the peri-apical infections of #1.2, 1.1 and 2.2.
5. Improve gingival status of inflammation.
6. Upgrade the quality of the restorations.
7. Keep it as conservative as possible.

TREATMENT ALTERNATIVES

Option 1: Conservative treatment plan with extraction of # 2.1

Eliminate old restorations in # 1.2, 1.1, 2.1 and 2.2

Retreat root canals treatments of # 1.2, 1.1 and 2.2

Ridge augmentation in #2.1 zone

Implant placement in #2.1 zone

Conventional fixed prostheses on # 1.2, 1.1 and 2.2

Screw retained implant-supported restoration in # 2.1

Replace old composite restorations in # 1.6, 1.7, 2.6, 2.7, 3.6, 3.7, 4.6 and 4.7

Option 2: Conservative treatment plan with extraction of # 1.1 and 2.1

Eliminate old restorations in # 1.2, 1.1, 2.1 and 2.2

Retreat root canals treatments of # 1.2 and 2.2

Ridge augmentation in #1.1 and 2.1 zone

Implant placement in #1.1 and 2.1 zone

Conventional fixed prostheses on # 1.2, and 2.2

Screw retained implant-supported restoration in # 1.1 and 2.1

Replace old composite restorations in # 1.6, 1.7, 2.6, 2.7, 3.6, 3.7, 4.6 and 4.7

PAGE 9

CLINICAL TREATMENT PLANNING

FINAL TREATMENT PLAN

One major goal driving the final plan was to keep it as conservative as possible. To that end, the following decisions were made:

Option 1: Conservative treatment plan with extraction of # 2.1.

Eliminate old restorations in # 1.2, 1.1, 2.1 and 2.2.

Retreat root canals treatments of # 1.2, 1.1 and 2.2.

Ridge augmentation in #2.1 zone.

Implant placement in #2.1 zone.

Conventional fixed prostheses on # 1.2, 1.1 and 2.2.

Screw retained implant-supported restoration in # 2.1.

Replace old composite restorations in # 1.6, 1.7, 2.6, 2.7, 3.6, 3.7, 4.6 and 4.7.

It was understood that this very conservative approach would result in several compromises:

There would be the possibility of failure of the root canals retreatments with its associated loss of the teeth involved. Also, #1.1 and its lack of facial bone are an unpredictable scenario of evolution, however the patient's age and requirements lead us to opt for this conservative treatment plan.

Initial therapy:

Initial therapy included oral hygiene instructions, periodontal prophylaxis and maintenance.

PAGE 10

ACTIVE CLINICAL TREATMENT



Pretreatment view showing # 1.2 full metal ceramic restoration. # 1.1, 2.1 and 2.2 relapsed composite veneers



Removal of the metal ceramic restoration, evidencing gingival inflammation and a non-precious alloy post and core.




Removal of composite veneers restoration, the image shows the minimum space between 2.1 and 2.2. This would be a critical issue in the future ridge augmentation, implant placement and restorative procedure.




Interim provisional restorations placed after the removal of relapsed restorations and before the surgical procedures.

PAGE 11

ACTIVE CLINICAL TREATMENT



Sequence of periapical radiographs showing the root canal retreatment of the maxillary incisor. It is highly recommended to interlace all the related disciplines to reach the best biological and functional possible outcome.



Post-endodontic intra-canal posts.

PAGE 12

ACTIVE CLINICAL TREATMENT



Surgical removal of the failing left central incisor. Note the absence of buccal plate due to the external resorption, and the lack of interproximal bone height between left and central incisors.



A large ridge defect is evident following the removal of the failing left central incisor. This scenario guide us to delay the implant placement and only proceed with the GBR with membrane barrier in this first surgical procedure.



Periodontal probe showing the absence of interproximal bone surrounding the mesial aspect of # 2.2



Horizontal and vertical ridge augmentation with 4BONE covered with a resorbable membrane.




PAGE 13


ACTIVE CLINICAL TREATMENT



Post surgical views 3 months following the GBR.



Post-surgical 3D ct scans after 7 months of graft healing



Pre-surgical 3D ct scans. Compare the volume successfully augmented

PAGE 14

ACTIVE CLINICAL TREATMENT



Set of photographs showing surgical procedure of the implant placement, implant selection MIS C1.

Successful orofacial ridge augmentation.

Good Interproximal bone in the distal aspect of the ridge.

Excellent interproximal bone in the mesial aspect of the ridge.

3-dimensional Ideal implant positioning.



One week post-surgical view.

Absence of distal peri-implant papillae

Lack of facial gingival tissue



PAGE 15

ACTIVE CLINICAL TREATMENT



Occlusal views of the esthetic zone showing the lack of facial gingival volume. The team decided to perform a new surgical procedure to gain soft tissue volume.





Soft tissue graft augmentation surgical procedure.



Clinical status 2 month after soft tissue surgery. Now we have the right scenario to proceed with the prosthetic restoration.

PAGE 16

ACTIVE CLINICAL TREATMENT



Set of photographs showing the switch from old composite restorations to new operative restorations in #1.6 and 1.7. This procedure appears to be simple and tedious, but we strongly believe in a full integration treatment plan, also it is well known that a correct approach in operative dentistry is the best path of preventive dentistry. This restorative procedures were made in the healing process of the soft tissue graft.

PAGE 17

ACTIVE CLINICAL TREATMENT



Set of photographs showing the switch from old composite restorations to new operative restorations in #2.6 and 2.7.

PAGE 18

ACTIVE CLINICAL TREATMENT



Set of photographs showing the switch from old composite restorations to new operative restorations in #3.6 and 3.7.

PAGE 19

ACTIVE CLINICAL TREATMENT



Set of photographs showing the switch from old composite restorations to new operative restorations in #4.6 and 4.7.

PAGE 20

ACTIVE CLINICAL TREATMENT

Provisionalization must be view as the non surgical refinement of the sof tissue architecture



New provisional restorations in place. Conventional fixed partial restorations in # 1.2 - 1.1 - 2.2.

Screw retained provisional restoration on implant in position of 2.1. The primary concern is to assure accurate fit of the crown margin to the implant shoulder, with no inclusion of cement, this is the reason why it is preferable to use an screw retained fixed restoration. Therefore it is critical to achieve the correct three dimensional position in the implant placement procedure as it result in the possibility of having the emergence axis of the screw in palatal position (cingulum). This costumized provisional restoration has the responsibility to shape the peri-implant gingival tissues. It is difficult to fully seat a definitive restoration if the peri-implant tissues have not been shaped with emergence profile provisional restorations.

As a result, there must be a precise transfer of this information to the technician about the clinician's ideal and the patient-approved soft tissue framework.

PAGE 21

ACTIVE CLINICAL TREATMENT



Before and after situation of peri-implant gingival tissue healing

PAGE 22

ACTIVE CLINICAL TREATMENT



Photograph showing clinical situation after 2 months of peri-implant gingival tissue management. Procedure of vital importance to ensure the right emergence profile of the future ceramic restoration.



Set of photographs detailing the simple and precise impression technique for the Mis C1 implant. It also shows the accurate impression of the prepared tooth. It is recommended for all ceramic restoration to finish the tooth preparation with a light chamfer.

PAGE 23

ACTIVE CLINICAL TREATMENT



Zirconium copyings try-in

You can also notice an anatomic wax-up in the implanted zone, used to certified and approve the emergence profile.



Set of photographs showing 4 All ceramic restorations.



Final All ceramic restoration in full harmony with the patients smile

ACTIVE CLINICAL TREATMENT



Before and after clinical situation 2 month installation.



PAGE 26



PAGE 27

CONCLUSION

From this clinical case discussion on implants ethics, it is clear that significant preplanning and an understanding of the various implant placement approaches and relative procedures have an important impact on minimizing negative hard and soft tissue contour changes from the moment of tooth extraction. One of the most important intent of this clinical case presentation is the intention to emphasize that implant treatment in the anterior maxilla zone requires a great knowledge base in both surgical and prosthetic aspects of treatment to meet esthetic ideal. It also underscores the key role of refinement of our treatment procedures, using more precise preoperative planning tools such as 3D computed tomography, three dimensional planning programs, ultra conservative but technically demanding surgical procedures, and absolute accuracy in provisionalization strategies. The result is that the previously well-defined roles of the surgeon and the restoring dentist are now more uncertain, emphasizing the benefits of developing implantologists or an implant team, that is, dentist or group of dentists who possess the knowledge and skills that are applied beginning planning phase, that demand excellence in the surgical phase, and that are full filled with artistry in the final restoration.

ACKNOWLEDGEMENTS

I would like to recognize the artistry of our dental ceramic technician Mr Angel Pricolo the enormous involvement of my implant and restorative team, Dr Hector Alvarez Cantoni and Dr Mariela Alvarez Cantoni.

Success and reaching for the 'patient's best' is thanks to our combined efforts.

PAGE 29

CURRICULUM VITAE OF THE AUTHOR

Joile Alvarez Cantoni

Graduate in year 2000 from Maimonides University, Buenos Aires, Argentina.

Associate professor for Partial Fixed prosthodontics department, Universidad de Buenos Aires, Argentina.

Chief of clinic for the Buccomaxillar rehabilitation postgraduate career, Universidad de Buenos Aires, Argentina.

Private dentistry in Cerrito 782 8 floor, Buenos Aires, Argentina (phone 54 11 4371 2971)

PAGE 30

REFERENCES

1. Buser D., Martin W, Belser UC. Optimizing esthetics for implant restorations in the anterior maxilla: Anatomic and surgical considerations. Int J Oral Maxillofac Implants 2004;19(supplement):43-61.
2. Schropp L, Wenzel A, Kostopoulos L, Karring T. Bone healing and soft tissue contour changes following single tooth extraction: A clinical and radiographic 12-month prospective study. Int J Periodontics Restorative Dent 2003;23:313-323.
3. Becker W, Goldstein M, Becker Be, Sennerby L. Minimally invasive flapless implant surgery: A prospective multicenter study. Clin Implant Dent Relat Res 2005;7 (supplement):S21-S27.
4. Bianchi AE, Sanfilippo F. Single tooth replacement by immediate implant and connective tissue graft: A 1-9 year clinical evaluation. Clin Oral Implants Res 2005;15:269-277.
5. Kinsel RP, Lamb RE. Tissue directed placement of dental implants in the esthetic zone for long-term biologic synergy: A clinical report. Int J Oral Maxillofac Implants 2005; 20:913-922
6. Tarrow D, Elian N, Fletcher P, et al. Vertical distance from the crest of bone to the height of the interproximal papilla between adjacent implants. J Periodontol 2003;74:1785-1788.
7. Chen ST, Wilson TJ Jr, Hammerle CHF. Immediate or early placement of implants following tooth extraction: Review of biologic basis, Clinical procedures, and outcomes. Int Oral J Maxillofac 2004;19(supplement):12-25.

PAGE 31

PATIENT'S & DOCTOR'S CONSENT

CONSENT TO ALLOW THE USE OF IMAGES AND OTHER RELEVANT TREATMENT DATA FOR PUBLICATION AND EDUCATIONAL PURPOSES

Dear Dr. Maclaren/Sir:

Please read the following information carefully. You are asked to sign this form as part of your submission of a case to the MIS Meeting Case Competition.

This form should be signed by both the treating dentist and by the relevant patient. Only forms including both signatures (D's and address) will be accepted.

If you have any questions, please contact Ms. Michael Molloy at MIS Implant Technologies LTD. michael@mis-implants.com

Please scan or send the signed form to michael@mis-implants.com, and please keep a copy of the form for your records.

PURPOSE

MIS Implant Technologies LTD is a company who manufactures and sells dental implants and other related components all over the world. As part of MIS's mission, the company aims to educate dentists and patients as for the benefits of dental implants. The company does in different ways to publish its products and to educate dental professionals as for its products and as for ways to provide best treatments to patients. In an attempt to enlarge our documented case library, MIS conducts a world-wide case competition, which will take place during MIS's Second Global Meeting in June 2013.

WHAT ARE YOU ASKED FOR?

Dentists: You are being asked to allow MIS Implant Technologies LTD to use images and other relevant treatment data that you submit for publication and for other activities, aiming to educate dental professionals and the public as for treatment done with the use of dental implants. The data that may be published or be used, in addition to the submitted images and treatment data will include specific recognition that it was provided by you. This recognition will be done by specifying your name, country and relevant data.

Patients: You are being asked to allow MIS Implant Technologies LTD to use images and other relevant treatment data for publications and for other activities, aiming to educate dental professionals and the public as for treatment done with the use of dental implants. The data that may be published or be used, in addition to the submitted images and treatment data will include: but not limited to your initials, age and gender, relevant data related to your health condition (especially if you are healthy or not, and if you suffer any conditions that may affect the success and outcomes of implant related treatments), smoking history, dental history, treatment plan, provided treatment and processes and outcome of treatments and more.

Although these images will be used without identifying information such as your name (other than your initials), and although measures will be taken to reduce or eliminate identifying features, the possibility remains that someone may recognize you.

PARTICIPANT STATEMENT

I (Dentist - Name and Identification number) JOSE P. MOLLOY ACU NO. 2620080
Address: CALLE 10 782 PISA, have read this consent form, and by signing this document I declare that:

1. I had an opportunity to ask questions, and if I did, I received answers that satisfied me completely.
2. I understand and that I allow MIS Implant Technologies LTD to use images, certain identifying information, and all other relevant submitted treatment data, for publication and for other activities, aiming to educate dental professionals and the public as for treatment done with the use of dental implants.
3. By signing this consent, I understand that I will not receive any direct or indirect compensation from MIS or any one of its representatives, in connection with this consent.

2-6-13

Date

Signature

I (Patient - Name and Identification number) ANITA JOSEFA MUAREZ DI: 75 361 950
Address: HOUBERTS 3224 2A, have read this consent form, and by signing this document I declare that:

1. I had an opportunity to ask questions, and if I did, I received answers that satisfied me completely.
2. I understand and that I allow MIS Implant Technologies LTD to use images, certain identifying information, medical information and all other relevant submitted treatment data, for publication and for other activities, aiming to educate dental professionals and the public as for treatment done with the use of dental implants.
3. By signing this consent, I understand that I will not receive any direct or indirect compensation from MIS or any one of its representatives, in connection with this consent.

2-6-13

Date

Signature

PAGE 32

CASE 3

EMILIO MATEO DOMINICANA



MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION
IMMEDIATE IMPLANT C1 MIS
WITH 4BONE GRAFT AND
IMMEDIATE LOADING IN THE
AESTHETIC ZONE.

Dr. Emilio Mateo

Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana

Consent and Release

mis

CONSENT TO ALLOW THE USE OF IMAGES AND OTHER RELEVANT TREATMENT DATA FOR
PUBLICATION AND EDUCATIONAL PURPOSES

Dear Dr./Madam/Sir:

Please read the following information carefully. You are asked to sign this form as part of your submission of a case to the MIS Meeting Case Competition.

This form should be signed by both the treating dentist and by the relevant patient. Only forms including both signatures ID's and address will be accepted.

If you have any questions, please contact Mr. Michal Malka at MIS Implant Technologies LTD: michal@mis-implants.com

Please scan or send the signed form to: michal@mis-implants.com, and please keep a copy of the form for your records.

PURPOSE

MIS Implant Technologies LTD is a company who manufactures and sells dental implants and other related components all over the world. As part of MIS's mission, the company aims to educate dentists and patients as for the benefits of dental implants. The company acts in different ways to publish its products and to educate dental professionals as for its products and as for ways to provide best treatments to patients.

In an attempt to enlarge our documented case library, MIS conducts a world-wide case competition, which will take place during MIS's Second Global Meeting in June 2013.

WHAT ARE YOU ASKED FOR?

Dentists: You are being asked to allow MIS Implant Technologies LTD to use images and other relevant treatment data that you submit, for publication and for other activities, aiming to educate dental professionals and the public as for treatment done with the use of dental implants. The data that may be published or be used, in addition to the submitted images and treatment data will include specific recognition that it was provided by you. This recognition will be done by specifying your name, country and relevant date.

Patients: You are being asked to allow MIS Implant Technologies LTD to use images and other relevant treatment data for publication and for other activities, aiming to educate dental professionals and the public as for treatment done with the use of dental implants. The data that may be published or be used, in addition to the submitted images and treatment data will include (but not limited to) your initials, age and gender, relevant data related to your health condition (especially if you are healthy or not, and if you suffer any conditions that may alter the success and outcomes of implant related treatments), smoking history, dental history, treatment plan, provided treatment and processes and outcome of treatments and more.

Although these images will be used without identifying information such as your name (other than your initials), and although measures will be taken to reduce or eliminate identifying features, the possibility remains that someone may recognize you.

PARTICIPANT STATEMENT

I (Dentist - Name and Identification number:) Emilio Mateo - 031-0967533-9

Address: Parque Mirador del Lago C/2 #17, have read this consent form, and by signing this document I declare that:

1. I had an opportunity to ask questions, and if I did, I received answers that satisfied me completely.
2. I understand and that I allow MIS Implant Technologies LTD to use images, certain identifying information, and all other relevant submitted treatment data, for publication and for other activities, aiming to educate dental professionals and the public as for treatment done with the use of dental implants.
3. By signing this consent, I understand that I will not receive any direct or indirect compensation from MIS or any one of its representatives, in connection with this consent.

15/11/12
Date

Dr. Emilio Mateo
Signature

I (Patient - Name and Identification number:) Wendy M. Baccaro 001-1019391-9

Address: Sancti Spiritus Hospital, Havana, have read this consent form, and by signing this document I declare that:

1. I had an opportunity to ask questions, and if I did, I received answers that satisfied me completely.
2. I understand and that I allow MIS Implant Technologies LTD to use images, certain identifying information, medical information and all other relevant submitted treatment data, for publication and for other activities, aiming to educate dental professionals and the public as for treatment done with the use of dental implants.
3. By signing this consent, I understand that I will not receive any direct or indirect compensation from MIS or any one of its representatives, in connection with this consent.

15/11/12
Date

Wendy M. Baccaro
Signature

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana


Preliminary Data

mis

A FULL SET OF PRELIMINARY RADIOGRAPHS



Initial
periapical
RX.



Sagittal CT
Scan

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Aravena
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana

MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT C.J.M.S WITH ABONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



mis

PRELIMINARY PHOTOGRAPHS



Intraoral and
extraoral
photographs

mis

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH ALBONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Astarup
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana

DETAILED SUMMARY OF DIAGNOSES

The clinical case corresponds to a 41-year-old patient, nonsmoker with no medical history of interest, who goes to the dental office by reference. After the clinical and radiological examination, it was diagnosed an agenesis of tooth # 23 with the presence of a composite restored primary canine. The patient reported that the restoration is often fractured, reason why she decided to find a most permanent solution. It was suggested upon her request an atraumatic extraction and immediate loading implant placement in the post-extraction socket.

TREATMENT OPTIONS

The options were:

- 3-piece fixed bridge.
- Two phases Implant.

-Post-extraction Implant and Immediate Loading:

The last option was the chosen one the Post-extraction Implant and Immediate Loading, because the patient rejected the 3 piece fixed bridge because she did not want "to sacrifice two healthy teeth" also she did not want two surgeries for the two phases implant.

DETAILED TREATMENT PLAN

After an effective treatment plan and computed tomography evaluation, it was place an infiltrative anesthesia in the area, the atraumatic extraction was done using periostomes and the integrity of the vestibular and palatal cortical was tested, following the principles of minimally invasive surgery, avoiding to lift flap or "flapless surgery", in order to preserve the soft tissue and minimize the resorption bone process.

The MIS C1 3.75 mm x 11.5 mm implant was placed, due to its conical anti-rotational connection with six positions and a position indicator, as well as its adaptation to abutment and the excellent seal that reduces micro-movements. The implant protocol was follow and the milling is ajusted to a more palatal axis instead of following the socket axis, placing the implant in a more palatalized position, which favors the regenerative process and preserve the buccal wall of the socket. The conical connection implant with platform switch, intended by design, is placed 2mm subcrestal. It is more efficient preserving the crestal bone and having a better impact on the formation of the papilla, which favors the anterior aesthetic.

Once the milling is finished, was proceed to the implant insertion using the contra-angle handpiece at low speed, until the implant was place in the ideal position in the three directions space, obtaining a primary stability over 60 N, which is an indispensable condition for the realization of an immediate restoration as planned in this case.

According to the current guidelines of the scientific evidence, we proceed to fill the GAP formed between the implant surface and the buccal wall of the socket, by the use of 4BONE, graft composed of 100% synthetic material, similar mineral human bone structure. The use of 60% hydroxyapatite with slow resorption rate and 40% of beta tricalcium phosphate with rapid resorption rate, guarantees a bone cell response in a manner similar to that caused by the bone, causing perfect balance. The 4BONE becomes a living bone with a new vascularization, due to angiogenic and osteogenic properties. This is introduced into the space using a syringe, taking care that the graft is not introduced in the implant interior using a healing screw.

For the temporary crown elaboration it was screwed a type PIK temporary abutment. From the confection of the diagnostic wax model, an acetate plate is made which was perforated through the acrylic provisional, it was adjusted without any contact and the acrylic crown made by the technician was filled with the autopolymerized acrylic.

Then the acetate matrix is removed where the provisional goes and the emerging profile of the crown was finishing and polishing, which is made in compliance with the dental organ dimensions, but eliminating the incisal edge to adjust occlusal surface and excursive movements to not induce micro-movements that could derail the implant osseointegration. This confirms that restoration contours that are in direct contact with the tissues must be fully polished, being placed the provisional prosthesis to the patient on the day of surgery.

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana

Treatment

mis

PHOTOGRAPHS



Preoperative image



Atraumatic extraction

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Aravena
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana

MIS



It verify if there are fenestration or dehiscence of the walls socket with a periodontal probe and locate the buccal ledge of the bone crest.



Milling the socket, flapless



Parallelism PIN



Placing the implant




Implant placed in the socket in a slightly palatalized position



MIS C1, 3.75mm. x 11.5 mm.
Implant

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Astarup
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Placement of immediate
temporary PIK



MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH ABONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Araven
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Confection of acrylic
temporary crown by the
prosthodontist from diagnostic
wax



MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH 4BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Aravena
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana




Inserting the healing screw
for bone graft placement
MIS 4BONE

MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH ABONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Aravena
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Remove the
healing screw for
placement of the
temporary crown.



MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH ABONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Aravena
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Temporary crown
p l a c e d
immediately after
surgery, once
removed the
cement and the
o c c l u s i o n
adjusted.

MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Immediate
post-surgical



One day after
post-surgical



7 days after
post-surgical

mis

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH ALONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana

Final Outcome

mis

SUMMARY OF THE CASE

After waiting 6 months for a good osseointegration process and after obtaining the main goals of an immediately post-extraction implant surgery with immediately loaded as: primary stability more than 60 N, the perfect balance to fill the gap with 4BONE graft and excellent adaptation of the gingival margin with the provisionally crown since the date of placement and without presenting any difficulty within the osseointegration period, we proceed to take the impression, with the open tray technique, often used in fixed prosthetic dentistry, with the double material technique with heavy and low to make the final crown.

The color was taken with the VITA Easyshade because it allows a digital reading that facilitates the clinical and laboratory work due to the evolution of dental colorimetric technology from which measurements are obtained color fast, accurate and objective.

We made all procedures required, and then we proceed to the placement of the final crown with a structure of CAD-CAM 3M Lava zirconia and ceramic VITA VM 9. Also, the necessary occlusal adjustments were made and the cement remains were removed to prevent inflammatory reactions in peri-implant tissues. Thus ends the post-surgical treatment, establishing a clinical and radiological follow up regularly to ensure the success of the implant.

SUMMARY OF THE PROPOSED CLINICAL PROTOCOL

Advances in implantology have grown as fast as the technology itself, reflected in planning atraumatic dental treatment and the use of precise guidelines that have led to the development of new techniques, allowing the insertion of dental implants and post-extraction minimally invasive surgery without the need of lifting a flap, with predictable results in function, aesthetics and patient comfort, keeping the original shape and thickness of the peri-implant tissues so that the post-surgical trauma and discomfort is minimizing compared to flap surgery, as long as these treatments accomplish the conditions and parameters for its realization.

That is why we propose a fully documented clinical case of a whole surgery with conical connection immediate implant and platform switch, the use of 4bone bone graft and immediately loaded in the esthetic zone, illustrating the alveolar bone preservation during and after extraction, which validates the atraumatic treatment advances and given evidence that this treatment is predictable.

RADIOGRAPHIC SEQUENCE



Initial Rx.



Parallelism PIN



Implant 2mm. subcrestal



Implant with temporal PIK

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Aravena
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana

MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Aravena
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Implant prosthetic
attachment



Crown final




Crown final.
Higher contrast radiography

MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Aravena
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana




Axial CT scan
Post - surgical



Sagittal CT scan
Post - surgical

mis

FINAL PHOTOGRAPHS



After 6 months the provisional is removed for making the final crown



Place the impression coping for open tray




Place resin flow to
copy the profile
emerging



Polymerizing the
resin flow



Open tray
technique




Printing technique of double material
heavy and silicone fluid



Replacing the analog

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Astarup
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Attachment prosthetic

MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH ABONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



It takes a
new interim
to improve
the
emergence
profile of
the canine
area



Profile
emerging
after the
placement of
the new
provisional

mis

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH ABONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Araven
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Structure in zirconium
3M LAVA CAD-CAM

mis

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Aravena
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana




Structure in zirconium

mis

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Printing drag

MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT C.J.M.S WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Astarup
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana




Structured final
zirconium crown 3M
LAVA CAD-CAM
and ceramic VITA
VM 9

mis

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (C.I.) MIS WITH ABONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Torque the final
attachment to 20 N,
is placed Teflon and
curing cement

MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Intraoral photos the
day of cementation

MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana



Intraoral photos the day of cementation

MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Internacional de la Rioja



Intraoral photos the day of cementation

MIS

Final photographs one month after of final restoration



Extraoral Photographs

QUIRURGIA INMEDIATA INVASIVA EN IMPLANTE INMEDIATO
POST-EXTRACCIÓN MISCL CON INIENTO 4BONE
Y CARGA INMEDIATA EN LA ZONA ESTÉTICA

Dr. Emilio Muñoz
Especialista en Periodoncia e Implantología Oral
Universidad Iberoamericana



MIS

QUIRURGIA INMEDIATA INVASIVA EN IMPLANTE INMEDIATO
POST-EXTRACCIÓN MISCL CON INIENTO 4BONE
Y CARGA INMEDIATA EN LA ZONA ESTÉTICA

Dr. Emilio Muñoz
Especialista en Periodoncia e Implantología Oral
Universidad Iberoamericana



mis



Intraoral Photographs

QUIRURGIA INMEDIATA INVASIVA EN IMPLANTE INMEDIATO
POST-EXTRACCIÓN MISCL CON INTENTO ABIGNE
Y CARGA INMEDIATA EN LA ZONA ESTÉTICA

Dr. Emilio Alvarado

Especialista en Periodoncia e Implantología Oral
Universidad Interamericana



MIS

QUIRURGIA INMEDIATA INVASIVA EN IMPLANTE INMEDIATO
POST-EXTRACCIÓN MISCL CON INLENTO 4BONE
Y CARGA INMEDIATA EN LA ZONA ESTÉTICA

Dr. Emilio Alonso

Especialista en Periodoncia e Implantología Oral
Universidad Iberoamericana



MIS

QUIRURGIA INMEDIATA INVASIVA EN IMPLANTE INMEDIATO
POST-EXTRACCIÓN MISCL CON INIENTO 4BONE
Y CARGA INMEDIATA EN LA ZONA ESTÉTICA

Dr. Emilio Alvarado
Especialista en Periodoncia e Implantología Oral
Universidad Iberoamericana



MIS

MINIMALLY INVASIVE SURGERY
FOR POST-EXTRACTION IMMEDIATE IMPLANT (CJ) MIS WITH BONE
GRAFT AND IMMEDIATE LOADING IN THE AESTHETIC ZONE.

Dr. Emilio Arias
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana

Acknowledgements

mis

Dr. Emilio Mateo (Periodontist and Implantologist case)

Dentist, Pontificia Universidad Católica Madre y Maestra.
Periodontics and Implant Dentistry, Universidad Iberoamericana.

Dr. Dioracy Vicioso (Rehabilitator of the case)

Dentist, Universidad Iberoamericana
Periodontics and Implant Dentistry, University of Sao Paulo, Brazil.

Elvin Santos.

Dental Technician, Universidad Iberoamericana.
SOLDESA Dental Laboratory.



Dr. Emilio Mateo
Specialist in Periodontics and Implant Dentistry
Universidad Iberoamericana

CASE 4

JAMES COLLINS DOMINICANA

**Implant Placement in the Esthetic Zone.
Surgical and Prosthetic Management.
Clinical Case Report**

Dr. James R. Collins C

Dr. Rubén T. Polanco A





1- Before extaction (frontal aspect of teeth 8 and 9)



2- Acrylic temporary restoration




3- Initial x-rays




4- 3D planification






5- The planned position for the future emergence are marked in the model



6- Platform switch titanium abutments are used as provisional



7- The titanium abutments are covered by resin




8- Peri implant area are also copy by composite




9- Acrylic temporary restoration finish before the surgery



10- Periodontal sounding to confirm the level of bone crest around adjacent teeth




11- Extraction of teeth 8 using a periotome instrument and a 150 forcep




12- Extraction of teeth 9 using a periotome instrument and a 150 forcep



13- Granulation tissue of the extraction socket is removed.



14. The blue circle shows the apical direction of the roots and the red is where the drilling will be done



15. The palatine inclination of the parallelism pins becomes clearly visible



17. Immediate temporal implant restoration on 8 and 9.



16. A sulcular incision is made with a no.15c blade



18. 1 month after surgery, clinical and radiographical evaluation. Note the use of platform switch abutments





19. 2 month after surgery, clinical and radiographical evaluation reveal stable periimplant tissue condition





20. During the temporal restorative phase, the development of aesthetic peri-implant soft tissue contours is achieved



21. Impression coping for the master cast.



22. Master cast.



23. Esthetic Zirconia abutments on 8 and 9.




24. Final abutments



25. A frontal view with the suprastructure in place for final details



26. An aesthetically pleasing overall integration of the two anterior reconstructions becomes clearly visible after its final cementation.



27. Occlusal view of the final restoration on 8 and 9 in the anterior maxilla



28. Before and 1 year clinical follow up




29. An aesthetically pleasing overall integration of the two anterior reconstructions is underlined by a close up view of the patient's non-forced smile.



CASE 5

M. L. RAMOS OLTRAS
SPAIN

LATERAL INCISORS AGENESIS



MARÍA L. RAMOS OLTRA

DDS, MSc, PhD

Introduction

- Tooth agenesis is one of the most common developmental dental anomalies.¹ Patients with congenitally missing teeth may present with undeveloped alveolar bone morphology, making implant reconstruction a challenge.² When the space maxillary isn't enough, treatment consisted of initial orthodontic space management to obtain adequate space for missing lateral incisors.³
- Dental Implants have become a primary treatment option for replacement of these teeth. Many times in prosthodontic treatment planning a multidisciplinary approach is needed for a comprehensive outcome. Prosthodontic treatment planning is needed prior to the patient's consultation and following treatment acceptance; the prosthodontist may need to coordinate treatment needs with other specialists, including an orthodontist and an implant surgeon.¹
- The demand for optimal orthodontic and prosthetic treatment is high because the condition has an impact on facial aesthetic.⁴

¹ Krossha M, Fickl S. Congenitally missing lateral incisors. A comparison between restorative, implant, and orthodontic approaches. Dent Clin North Am. 2011 Apr;55(2):283-99.

² Nilsson J, Mordinaer O, Strauss M, Pelea M, Sacco P, Chausu G. Implant-supported restoration of congenitally missing teeth using cancellous bone block-allografts. Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2011 Mar;111(3):286-91.

³ Patil PG, Koremore V, Chavan S, Nimbalkar-Patil SB, Kulkarni B. Multidisciplinary treatment approach with one piece implants for congenitally missing maxillary lateral incisors: a case report. Eur J Prosthodont Restor Dent. 2012 Jun;20(2):92-6.

⁴ Robertsson S, Mohlin B, Thilander B. Aesthetic evaluation in subjects treated due to congenitally missing maxillary laterals. A comparison of perception in patients, parents and dentists. Swed Dent J. 2010;34(4):177-86.

Treatment planning

- ✦ This case describes multidisciplinary management presenting thin spaced maxillary anteriors due to the congenitally missing lateral incisors.
- ✦ Single piece in 12 and 22 position, standard diameter implants (3,75 x 13mm C1 MIS®) were placed in edentulous spaces on both sides. During the surgery, ROG was made using MIS® 4BONE (Synthetic Bone Graft 0,5cc, particula size 0,5-1mm) and Connective tissue graft in 22 position.
- ✦ Resin crowns were given as provisional restorations.
- ✦ Metal-ceramic crowns were given as definitive restorations, resulting into an acceptable aesthetic outcome.

Conclusion

- In the rehabilitation of a single missing lateral maxillary incisor, no statistically significant difference was assessed between immediately and one-stage restored small-diameter implants with regard to implant survival, mean marginal bone loss, and probing depth. Narrow implants proved to be a predictable treatment option if a strict clinical protocol was followed.⁵
- Successful and satisfying dental treatment is always the goal for patients and dental practitioners, meaning that a patient's needs are solved in a functional and esthetic way. Patients and dentists have to find the best way to reach their common goal of satisfaction. Some authors introduces examples of different approaches to solve the problem of congenitally missing lateral incisors. In most cases, an interdisciplinary treatment plan has to be worked out and executed.¹

⁵ DeGidi M, Nardi D, Pinelli A. Immediate versus one-stage restoration of small-diameter implants for a single missing maxillary lateral incisor: a 3-year randomized clinical trial. J Periodontol. 2009 Sep;80(9):1393-8.

¹ Krassig M, Eickl S. Congenitally missing lateral incisors. A comparison between restorative, implant, and orthodontic approaches. Dent Clin North Am. 2011 Apr;55(2):283-99.

Intraoral Analysis



Sex: Female
Date of birth: 1992

June 2012

Intraoral Analysis



June 2012

Intraoral Analysis



June 2012

Intraoral Analysis



Rx



CT scan




Risk Factors


	DANGER	CAUTION	OKAY
GENERAL HEALTH			
			Good
PATIENT INTERROGATION			
Esthetic demands			
Motivation			
ETIOLOGY OF THE EDENTULISM			
			Agensis
EXTRAORAL EXAMINATION			
Smile line		Medium	

INTRAORAL EXAMINATION			
Jaw opening			Good
Hygiene			Good
Lesions, abscess			No
Vestibular concavity present			No
Discrepancy mx-mdb.			No
Intraoral palpation			Good
Vertical bone resorption			
Height between bone crest and opposing tooth			
Gum			
Papillae of adjacent teeth			
Dental risks factors			
Provisionalization	Immediate		

FUNCTIONAL EVALUATION			
Bruxism/parafunction			No
Natural teeth participating in proprioception			Yes
Lateral contacts			
RADIOGRAPHIC EXAMINATION			
Chronic lesions close to the implant zone			No
PERIODONTAL EVALUATION			
Gingivitis			No
Treated periodontitis			No
Active periodontitis			No



Surgical Procedure 12i






June 2012

Surgical Procedure 22i





✓ Connective tissue graft



June 2012


12





22



CT scan



4Bone



Provisionalization



July 2012

Final Crowns



Acknowledgements:

- Dental technician PACO CARRILLO (CLASSIC DENTAL, Murcia, SPAIN).
 - Prof. JOSE LUIS CALVO GUIRADO. (Faculty of Dentistry, Murcia University, Murcia, SPAIN).
-

CASE 6

M. P.

RAMÍREZ FERNÁNDEZ
SPAIN



Concurso de Casos Clínicos Mis II Congreso Cannes Francia 2013



Dra. María Piedad Ramírez Fernández Dds, Msc, Phd,

Universidad De Murcia

mpramirezfern@hotmail.com

Teléfono: 630634344

Datos preliminares

Historia clínica

DATOS DE FILIACIÓN

Nombre: Joaquín Templado Gómez

Sexo: Varón

Edad: 59

Ocupación: Profesor

Domicilio: C/ Mesones, 10, 2º C (Cieza) Murcia

Estado Civil: Casado

ANTECEDENTES PERSONALES:

Intervenciones Quirúrgicas: Reducción de estómago

Alergias: No Conocidas

Hábitos: Fumador 20 cigarrillos / día

ANTECEDENTES FAMILIARES:

Familia Materna: Diabetes y Cáncer e Infarto de Miocardio.

ANAMNESIS

MOTIVO PRINCIPAL DE CONSULTA: Restaurar la dentición con implantes dentales.

EXÁMENES COMPLEMENTARIOS

-I cat vision

DIAGNÓSTICO

Paciente parcialmente edéntulo con dientes naturales remanentes a nivel anterior.

Caries múltiple a nivel radicular: 12, 11, 21, 22, 35, 34, 33, 43, 44.

Movilidad dentaria a nivel: 31, 32, 41, 42.

Múltiples restos radiculares: 14, 13, 23, 24, 25, 26, 27, 45

Ausencia 15, 16, 17, 18, 36, 37.

Quiste residual a nivel del 36, 37.

PLAN DE TRATAMIENTO

Plan de tratamiento A

-Exodoncias completa de piezas dentarias de arcada superior e inferior.

-Extirpación de lesión quística y regeneración ósea del defecto.

-Colocación de 8 implantes en el maxilar superior y 8 en el maxilar inferior para realizar tratamiento con prótesis fija cementada sobre implantes.

Plan de tratamiento B

-Exodoncias completa de piezas dentarias de arcada superior e inferior.

-Extirpación de lesión quística y regeneración ósea del defecto.

-Colocación de 6 implantes en maxilar superior (All on six) y 4 implantes en mandíbula (All on four) para una prótesis híbrida atornillada.

Plan de tratamiento C

-Exodoncias completa de piezas dentarias de arcada superior e inferior.

-Extirpación de lesión quística y regeneración ósea del defecto.

-Colocación de 4 implantes en maxilar superior y 2 implantes en mandíbula para realizar una rehabilitación con sobredentaduras.

PLAN DE TRATAMIENTNO ELEGIDO

Considerando la edad del paciente y las características del hueso remanente elegimos el plan de tratamiento A, el cual le permitía al paciente una mayor calidad de vida.

Tratamiento

1. Exodoncia completa de piezas dentales remanentes y restos radiculares.
2. Extirpación de quiste residual y regeneración ósea utilizando una mezcla de hueso Bond Bone y 4Bone para regenerar el defecto.
3. Colocación de Implantes dentales inmediatos post-extracción.
MAXILAR SUPERIOR: 11,13, 15, 16, 21, 23, 25, 26 IMPLANTES C1- MIS
MANDIBULAR: - 31, 33, 35, 36 IMPLANTES SEVEN-MIS
- 41, 43, 45, 46 IMPLANTES C1- MIS
4. Colocación de 6 micro-implantes transicionales en maxilar superior y 4 micro-implantes transicionales en mandíbula para soportar los provisionales
5. Configuración de juego de provisionales en resina.
6. Radiografía de control
7. Pasados 5 meses de evolución, levantamos los provisionales
8. Retiramos los micro-implantes transicionales a excepción de los dos más posteriores del maxilar que servirán a la provisionalización.
9. Descubrimiento de los implantes sumergidos.
10. Colocación de transfer de impresión
11. Ferulización con alambre de ortodoncia de los pilares
12. Sellado con resina autopolimerizable
13. Toma de impresión a cubeta abierta con silicona.
14. Colocación de pilares provisionales para provisionalización.
15. Registro de las relaciones craneomaxilares y toma de la dimensión vertical.
16. Prueba de los pilares en boca y ferulización con resina autopolimerizable
17. Colocación de los pilares definitivos.
18. Cementado de la estructura fija metal-porcelana.
19. Radiografía de control al mes de la colocación de prótesis definitiva



Ilustración 1. Imagen preoperatoria del paciente.



Ilustración 2. Radiología preoperatoria.



Ilustración 3. Exodoncias múltiples



Ilustración 4. Implantes C1- MIS



Ilustración 5. Implantes Seven-Mis




Ilustración 6. Colocación de 4 implantes seven lado izquierdo mandibular.




Ilustración 7. Colocación de 8 implantes C1 a nivel maxilar y 4 implantes C1 en la do derecho mandibular.




Ilustración 8 . Extirpación de Quiste Residual a nivel mandibular.



Ilustración 9. Bond Bone



Ilustración 10. 4 Bone



Ilustración 11. Mezcla de Bond- Bone y 4 Bone para regenerar el defecto óseo.



Ilustración 12. Preparación del rehidratado del 4 Bone



Ilustración 13. Relleno del defecto óseo y del gap en implantes inmediatos con una mezcla de Bond-Bone y 4 Bone.



Ilustración 14. Cierre completo de la cirugía con suturas dobles y simples.



Ilustración 15. Radiografía de control de la cirugía.



Ilustración 16. Provisionalización inmediata.



Ilustración 17. Imagen postoperatoria 5 meses de evolución.



Ilustración 18. A bordaje para configuración de la prótesis definitiva.



Ilustración 19. Implantes definitivos sumergidos.



Ilustración 20. Retirada de microimplantes a nivel maxilar.




Ilustración 21. Retirada de microimplante en mandíbula



Ilustración 22. Descubrimiento de los implantes sumergidos a nivel del maxilar. Colocación de pines de transferencia implante C1 en maxilar.



Ilustración 23. Ferulización de pilares



Ilustración 24. Ferulización con resina




Ilustración 25. Toma de impresiones maxilar.



Ilustración 26. Colocación de transfer a nivel del maxilar inferior



Ilustración 27. Ferulización con resina



Ilustración 28. Pilares provisionales



Ilustración 29. Adaptación de los provisionales.



Ilustración 30. Control radiológico.



Ilustración 31. Toma de la dimensión vertical



Ilustración 32. Prueba de pilares definitivos pasados 15 días de la toma de impresiones.



Ilustración 33. Ferulización de pilares definitivos.



Ilustración 34. Llaves de transferencia para la colocación de los pilares definitivos.



Ilustración 35. Llave dinamométrica para realizar el atornillado definitivo de los pilares.



Ilustración 36. Colocación de pilares definitivos pasados 10 días de la prueba de pilares a nivel maxilar.



Ilustración 37. Colocación de pilares definitivos y cierre de las chimeneas pasados 10 días de la prueba de pilares a nivel mandibular.



Ilustración 38. Vista oclusal de la estructura metal-porcelana a nivel mandibular.



Ilustración 39. Vista dorsal de la estructura previa al cementado

4. Resultado Final



Ilustración 40. Vista panorámica de los resultados .



Ilustración 41. Vista frontal a mayor aumento



Ilustración 42. Vista del frente estético a nivel frontal.



Ilustración 43. Vista de la zona mucogingival.



Ilustración 44. Vista lateral de la zona estética.



Ilustración 45. Resultado final del caso clínico



Ilustración 46. Radiología de seguimiento al mes de la colocación

5. Reconocimientos

Agradecemos la colaboración del taller de prótesis Clasic Dental en Murcia.

CASE 7

E. DUPERLY SANCHEZ COLOMBIA

Single tooth replacement and immediate loading at the esthetic zone with C1 MIS implant: a multidisciplinary approach.

Eduardo Duperly Sanchez D.D.S.*, Agnes Kristine Lanner D.D.S.**, Maria Lorena Yepes Sáenz D.D.S.***

Dental implants are a very good choice to replace a missing tooth, however sometimes there are not the optimal conditions to make a post extraction implantation. As a result, procedures like socket preservation or improvement and guided bone regeneration (GBR) become necessary.

A 38 year old male patient came complaining of pain in his tooth 9 (fig. 1), the clinical examination shows a palatal cervi-

cal inflammation. An increased probing depth of 12 mm on the palatal aspect was shown. At the periapical radiography was an apparent internal root resorption (fig. 2). He was referred to root canal where a communication between pulp space and the external portion of the root at the palatal aspect was evidenced.

To get a more accurate diagnosis he was asked to make a CT scan, where a communication to the palatal aspect and an apparent line of fracture at the buccal medial third of the root appeared (fig. 3); three alternatives of treatment were proposed:

1. Extraction, socket preservation and implantation 6 months later with immediate loading, and a complete fixed orthodontic treatment.

2. Extraction, socket preservation and implantation 6 months later without immediate loading and a complete fixed orthodontic treatment.

3. Extraction, immediate implantation and GBR without immediate loading and a complete

fixed orthodontic treatment.

The team decided the first treatment plan.

The patient was informed about all relevant aspects of the treatment. Based on this comprehensive information, he agreed the proposed treatment and gave a written informed consent.

The first treatment step was the careful extraction. This was carried out without flap elevation (fig. 4). The extraction socket was carefully debrided and filled with a 0.5 cc Puros Cancellous Particulate Allograft small particle (Zimmer®) and Socket Repair Membrane (Zimmer®, fig. 5). Once extracted, the tooth was sectioned horizontally at the crevical root level to get a closer look of the communication (fig. 6). In addition, The patient needed also an orthodontic treatment to improve occlusion and to make the closure of some dental spaces that were left from a previous treatment. After that, fixed orthodontic appliances were put with two purposes:

1. To hold a temporary or pro-

*D.D.S. Prosthodontist, Assistant Professor, Department of Prosthodontics, Faculty of Dentistry Javeriana University Bogotá D. C., Colombia

** D.D.S. Orthodontist, Private practice. Faculty of Dentistry Javeriana University Bogotá D. C., Colombia

***D.D.S. Periodontist, Periodontist, Instructor professor, Department of Periodontics Faculty of Dentistry Javeriana University Bogotá D. C., Colombia

visional, and 2. To solve occlusal and esthetic problems.

Without making pressure on the soft tissues, the provisional was put "hanging" between the adjacent teeth. After a healing period of 6 months the clinical examination revealed uneventful healing of the extraction socket, and adequate position of the teeth in terms of required space

for implantation and restoration. The soft tissues were fully intact (fig. 7).

Preoperative analysis of the anatomical conditions was evaluated using a CT scan, and the prosthetic planification of the surgery and surgical guide was carefully done (fig. 8, 9, 10).

The second procedure was performed passed 6 months with

a fullthickness flap raised and a conical connection MIS implant C1 (3.75 x 13mm) was placed confirming a final torque of 40 Ncm (Fig. 11, 12). Special emphasis was placed on obtaining a correct three-dimensional position with the help of the surgical guide. Apical implant threads were exposed on the buccal surface. Mis Bond Bone™ was used to cover



Fig. 1 Initial situation



Fig. 2
Initial
X ray




Fig. 3 Initial CT scan



Fig. 4 Probing after extraction



Fig. 5 Socket repair technique



Fig. 6 Cross section of tooth



Fig. 7 6 months after extraction




Fig. 8 CT scan pre implantation



Fig. 9 Wax up of tooth



Fig. 10 Surgical guide



Fig. 11 Drilling



Fig. 12 C1 MIS implant™



Fig. 13 Apical fenestration



Fig. 14 Bond Bone MIS™



Fig. 15 Bond Bone covering



Fig. 16 Peek post in position



Fig. 17 Composite temp. crown



Fig. 18 Repositioned flap



Fig. 19 Gingival contouring



Fig. 20 After 4 months of healing

the exposed threads of the implant, offering the osteoconductive capability required to achieve regeneration of vital bone (Fig. 13, 14, 15).

With the implant in adequate position, the temporary Peek post was used to make immediate provisionalization developing a smooth composite emergence profile of the incisor required to

enhance gingival esthetics (Fig. 16, 17). Repositionary of the flap was achieved and tension free utilizing 6-0 Prolene Ethicon suture (Fig. 18).

After 4 months, once the case achieved the complete healing and tissue maturation, (Fig. 19, 20) impression was taken with polyvinyl siloxane impression material. A titanium base abutment

was used to make a Prettau Zirconia cemented crown; final x ray control was taken. (Fig. 21, 22, 23).

The treatment of the case purposed from the beginning achieved its objectives, such as preservation of crestal bone and as a consequence, gingival esthetics.



Fig. 21 Final tissue profile



Fig. 22 Final restoration



Fig. 23 Final X ray control

REFERENCES

- 1.- Momen A. Atieh, Hadeel M. Ibrahim, and Ahmad H. Atieh "Platform Switching for Marginal Bone Preservation Around Dental Implants: A Systematic Review and Meta-Analysis". *J Periodontol* 2010;81:1350-1366
- 2.- Maria Bateli, DDS/Wael Att, DDS, Dr Med Dent Habil/Jörg R. Strub, DDS, Dr Med Dent, Dr hc, PhD. "Implant Neck Configurations for Preservation of Marginal Bone Level: A Systematic Review". *Int J Oral Maxillofac Implants* 2011; 26:290-303
- 3.- Nicolas Elian, DDS, Sang-Choon Cho, DDS, MS, Stuart Froum, DDS, Richard B. Smith, DDS, Dennis P. Tarnow, DDS. "Simplified Socket Classification and repair technique". *Practical Procedures & Aesthetic Dentistry* 2007; 19(2):99-104
- 4.- B. Shakibaie-M. "Socket and ridge preservation from the three-dimensional perspective". *Sonderdruck Zeitschrift für Zahnärztliche Implantologie* 25. Jahrgang / Heft 4/2009, Seiten 369-377
- 5.- Stephen L. Wheeler, DDS "Implant Complications in the Esthetic Zone". *J Oral Maxillofac Surg* 65:93-102, 2007, Suppl 1