

Internal Hexagon Prosthetic Options



mis[®] | MAKE IT SIMPLE

The MIS Quality System complies with international quality standard ISO 13485: 2016 - Quality Management System for Medical Devices, with Medical Device Directive 93/42/EEC, and with EU Medical Device Regulation MDR 2017/745.

IFUs for MIS products may be found at: <https://ifu.mis-implants.com>. Adobe Acrobat is required to view the IFU file on the website. This software may be freely downloaded from the Adobe website

Table of Contents

Implant Dimensions

P. 6 - 7

Prosthetic Options Listed by Healing Caps

P. 8 - 9

Healing Caps

P. 10 - 12

Impression Coping

P. 13 - 14

Temporary Restoration

P. 15 - 17

○ Temporary Cylinders

Screw-Retained Restoration

P. 18 - 23

- Multi-Unit System
- Plastic Cylinders

Digital Restoration

P. 24 - 29

- EZ-Base System
- Ti-Bases
- Digital Model Analogs
- Scan Posts
- Titanium Blanks

Cemented Restoration

P. 30 - 39

- Cementable Abutments
- CPK

Overdenture

P. 40 - 45

- LOCKiT System
- OT-Equator Kits
- Ball Attachment System

Prosthetic Screws and Drivers

P. 46 - 47

Implant Dimensions: SEVEN

Narrow
Platform







Standard
Platform



Standard
Platform



Wide
Platform




Implant Dimensions: M4




Narrow
Platform



Standard
Platform



Wide
Platform




ALL PLATFORMS

Prosthetic Options Listed by Healing Caps


Each column lists the different heights of healing caps and corresponding prosthetic options. **Bold** item numbers indicate the products shown in the images above.

*Image is for illustration only. (2mm height closed tray impression coping presented. Not available in 3mm).


1. Narrow Platform



2. Standard Platform



3. Wide Platform



1.



	Healing caps	Impression coping open tray	Impression coping closed tray	Cemented options	Digital options
	4mm	3.2mm	4mm	4mm	4.2mm
Tissue Level	2mm				
Bone Level	3mm				
	MH-N0440	MN-IO401	MN-ITC30*	MK-NPC62 MN-MC102 MN-AT102 MN-AT202	
	MH-N0540			MK-NPC63	MN-TBC30 MN-FTB30
	MH-N0640			MK-NPC64	
Gingival heights 2, 3, 4mm					

2.

Healing caps

Impression coping open tray

Impression coping closed tray

Cemented options


Screw retained options

Digital options

or

or

or



MH-S0448

MD-IO481

MD-ITC30*

MK-CPC62
MD-MC102
MD-AT152
MD-AT252

MA-SC248

MH-S0548

MD-AT153
MD-AT253
MK-CPC63

MA-SC348

MD-TBC30
MD-FTB30

MH-S0648

MK-CPC64

MA-SC448

Gingival heights 2, 3, 4mm

3.

Healing caps

Impression coping - open tray

Impression coping - closed tray

Cemented options


Screw retained options

Digital options

or

or

or



MH-W0455

MW-IO551

MW-ITC30*

MK-WPC62
MW-MC102
MW-AT152

MH-W0555

MK-WPC63

MA-WC348

MW-TBC30
MW-FTB30

MH-W0655

MK-WPC64

Gingival heights 2, 3, 4mm


Healing Caps Narrow Platform

i Tightening
torque:
20-25Ncm

Thin

Ø4

H2




MH-N2330

H3




MH-N3330

H4




MH-N4330

H5



MH-N6330

H6



MH-N8330

Standard


Ø4



MH-N0440




MH-N0540




MH-N0640

Anatomic


Ø4.3




MH-52330




MH-53330



MH-54330



MH-55330



MH-56330

Healing Caps Standard Platform

i
 Tightening torque:
 20-25Ncm

Thin

Ø4

Standard

Ø4.8

Standard

Ø4.8

CONCAVE

Anatomic

Ø5.5

H3

H4

H5

H6



MH-03375



MH-04375



MH-05375



MH-06375



MH-A3375




MH-A4375




MH-A5375



MH-A6375



MH-S0448



MH-S0548



MH-S0648



MH-53375



MH-54375



MH-55375



MH-56375

Healing Caps Wide Platform

i Tightening
torque:
20-25Ncm

Standard Ø5

Standard Ø5.5



Anatomic Ø6.5

H2

H3

H4

H5


H6



MH-W3500



MH-W4500



MH-W5500



MH-W0455



MH-W0555




MH-W0655



MH-W3630



MH-W4630



MH-W5630

Impression Copings Open tray

i Tightening torque: 20-25Ncm

Impression coping

Impression coping










CONCAVE

Implant analog,
open/closed tray

NARROW

STANDARD

WIDE



Impression Copings

Closed tray

i Tightening torque: 20-25Ncm

Impression coping

Impression coping set


CONCAVE

Press fit


NARROW

STANDARD


WIDE




MN-IT100




MD-IT100




MW-IT100




MN-IT300




MD-IT300




MW-IT300




MN-ITC30




MD-ITC30




MW-ITC30



MN-PF330



MD-PF375



MW-PF550

Temporary Restoration

Temporary Plastic Cylinders Anti-rotation

i
Tightening
torque:
20-25Ncm

Narrow profile temporary
plastic abutment

H1

Temporary plastic
abutment

H1


Temporary
plastic cylinder

H2.5

NARROW

STANDARD

WIDE




MN-TPA12




MD-TPA12




MW-TPA12




MN-TPA10




MD-TPA10



MW-TPA10



MD-TPH50



MW-TPH50

Temporary Restoration


Temporary Cylinders

Free rotation

i Tightening
torque:
20-25Ncm

H1

NARROW




MN-TC013

STANDARD



MD-TC013

WIDE




MW-TC013

H2

CONCAVE



MN-TC023



MD-TC023




MW-TC023

H3


CONCAVE



MN-TC033



MD-TC033



MW-TC033

Temporary Restoration

Temporary Cylinders Anti-rotation

i
Tightening
torque:
20-25Ncm

H1

H2

CONCAVE


H3

CONCAVE


NARROW

STANDARD


WIDE




MN-TCH13



MD-TCH13




MW-TCH13




MN-TCH23




MD-TCH23




MW-TCH23



MN-TCH33



MD-TCH33



MW-TCH33


Multi-Unit System

Tightening torque: 30Ncm




H2

NARROW



MA-NC248


STANDARD




MA-SC248

WIDE

H3



MA-NC348




MA-SC348




MA-WC348

H4




MA-NC448



MA-SC448

H5



MA-WC548

Screw Retained Restoration

Multi-Unit System

Tightening torque:
30Ncm

H1

NARROW



MA-N1480

STANDARD




MA-S1480

WIDE




MA-W1480

H2




MA-N2480




MA-S2480


H3



MA-N3480



MA-S3480




MA-W3480

H4




MA-N4480



MA-S4480


H5



MA-N5480



MA-S5480



MA-W5480

Screw Retained Restoration


Multi-Unit System

Angulated abutments,
17°

i
Tightening
torque:
20Ncm


H1

NARROW




MA-NO171

STANDARD




MA-SO171

WIDE




MA-WO171


H2



MA-NO172




MA-SO172




MA-WO172

H3



MA-SO173



MA-WO173

Screw Retained Restoration

Multi-Unit System


Angulated abutments,
30°

i

Tightening
torque:
20Ncm


H1

NARROW




MA-NO301

STANDARD




MA-SO301

WIDE




MA-WO301


H2



MA-NO302




MA-SO302




MA-WO302

H3




MA-SO303



MA-WO303


Multi-Unit System Components

Healing cap




MU-H4480

Analog




MU-RSM48

Impression coping,
open tray




MU-I0480

Impression coping,
closed tray




MU-IT480

Temporary cylinder



MU-TO480

Plastic cylinder




MU-CP048

Gold bar coping abutment



MU-BGC48

Gold plastic cylinder



MU-MG480

CAD/CAM cementing cap (New design)



MU-CC481

CAD/CAM cementing cap



MU-CC480

Scan Post




MU-SP102

Digital model analog




MU-MAN10

Straight Multi-Unit
ratchet key



MT-MURL2 - Long
MT-MURS2 - Short

Straight Multi-Unit
motor key



MT-MUML2 - Long
MT-MUMS2 - Short

Screw Retained Restoration

Plastic Cylinders

Free rotation

i Tightening
torque:
30Ncm


Gold plastic
cylinder

NARROW




MN-GP010

STANDARD




MD-GP010

WIDE




MW-GP010


Plastic cylinder



MN-CP013



MD-CP013




MW-CP013

CoCr plastic
cylinder



MN-CC010




MD-CC010




MW-CC010

Wide profile
plastic cylinder



MD-CP050



MW-CP060


Screw Retained Restoration

Plastic Cylinders Anti-rotation


Tightening
torque:
30Ncm

Gold plastic
cylinder


NARROW






STANDARD






WIDE





Plastic cylinder



CoCr plastic
cylinder



Wide profile
plastic cylinder



Digital Restoration

EZ-Base System




Tightening torque:
25Ncm




H1.5


NARROW



STANDARD



WIDE




EZ-Base System Components

Ratchet key for angled screw channel




MT-ELR10 - Long
MT-ESR10 - Short

Motor key for angled screw channel




MT-ELM10 - Long
MT-ESM10 - Short

EZ-Base prosthetic screw



MN-E0160 - NP
MD-E0220 - SP






A unique screw and driver tip design enable extreme angulation and a wide insertion path.




Ti-Base
Anti-rotation

i
Tightening
torque:
30Ncm




H0.5





H1.5



H3



H0.5



Ti-Base
Free rotation




Tightening
torque:
30Ncm

H0.5




NARROW

STANDARD



WIDE



H1.5



H3



Digital Restoration

Digital Model Analog

Tightening
torque:
15Ncm

NARROW




STANDARD



WIDE



MULTI-UNIT



Scan Posts



Tightening
torque:
20Ncm

NARROW

STANDARD

WIDE

MULTI-UNIT




MN-SP102



MD-SP102




MW-SP102



MU-SP102

MK-0056 SCAN POST KIT, INTERNAL HEX.



MIS scan posts are also sold in a kit.

Digital Restoration

Titanium Blanks Anti-rotation

i
Tightening
torque:
30Ncm

Ø9

Ø12

Amann Girrbach

Ø12

Milling analogs

NARROW

Ø9mm



MN-IB091

STANDARD

Ø9mm



MD-IB091

WIDE

Ø9mm



MW-IB091

Ø12mm



MN-IB121

Ø12mm



MD-IB121

Ø12mm



MW-IB121

Ø12mm



MN-AGT12

Ø12mm




MD-AGT12

Ø12mm




MW-AGT12

Ø5




MN-MAB10

Ø5



MD-MAB10

Ø5



MW-MAB10

Cemented Restoration


Cementable Abutments

i Tightening
torque:
30Ncm



0°

NARROW


STANDARD












10°



15°



WIDE





Cemented Restoration

Cementable Abutments

i Tightening
torque:
30Ncm

20°







25°


NARROW

STANDARD


Cementing posts




STANDARD




MD-A2510




MD-AN251



MD-P2530




MD-AT252




MD-AT253


WIDE




MD-CTP10




MD-MACF1



MD-MC102




MW-MAC10



MW-WMAC1



MW-CTP10





MW-MC102

Cemented Restoration

CPK - Transgingival Abutments

i Tightening
torque:
30Ncm







NARROW

STANDARD


WIDE




MN-CPC62




MD-CPC62




MW-CPC62




MN-CPC63




MD-CPC63




MW-CPC63



MN-CPC64




MD-CPC64




MW-CPC64

STANDARD


WIDE




MD-CPK81




MW-CPK81




MD-CPK82




MW-CPK82




MD-CPK83



MW-CPK83



MD-CPK84



MW-CPK84

CPK - Kit Contents

Burnout plastic cap,
anti-rotation




MN-ICH40
MD-ICH40
MW-ICH40

Burnout plastic cap,
free rotation




MN-IC040
MD-IC040
MW-IC040

Plastic healing cap




MM-CHC40
MM-CHC60
MM-CHC80

Transfer plastic cap



MN-IC600
MD-IC800
MW-IC800

Analogs



MN-RSM60
MD-RSM60
MW-RSM60

CPK abutments are
also sold in kits.



MK-NPK41 ● Ø4, H=1, C=4
MK-CPK41 ● Ø4.8, H=1, C=4
MK-WPK41 ● Ø5.5, H=1, C=4

MK-NPK42 ● Ø4, H=2, C=4
MK-CPK42 ● Ø4.8, H=2, C=4
MK-WPK42 ● Ø5.5, H=2, C=4

MK-NPK43 ● Ø4, H=3, C=4
MK-CPK43 ● Ø4.8, H=3, C=4
MK-WPK43 ● Ø5.5, H=3, C=4

MK-NPK44 ● Ø4, H=4, C=4
MK-CPK44 ● Ø4.8, H=4, C=4
MK-WPK44 ● Ø5.5, H=4, C=4

MK-NPK61 ● Ø4, H=1, C=6
MK-CPK61 ● Ø4.8, H=1, C=6
MK-WPK61 ● Ø5.5, H=1, C=6

MK-NPC62 ● Ø4, H=2, C=6
MK-CPC62 ● Ø4.8, H=2, C=6
MK-WPC62 ● Ø5.5, H=2, C=6

MK-NPK62 ● Ø4, H=2, C=6
MK-CPK62 ● Ø4.8, H=2, C=6
MK-WPK62 ● Ø5.5, H=2, C=6

MK-NPC63 ● Ø4, H=3, C=6
MK-CPC63 ● Ø4.8, H=3, C=6
MK-WPC63 ● Ø5.5, H=3, C=6

MK-NPK63 ● Ø4, H=3, C=6
MK-CPK63 ● Ø4.8, H=3, C=6
MK-WPK63 ● Ø5.5, H=3, C=6

MK-NPC64 ● Ø4, H=4, C=6
MK-CPC64 ● Ø4.8, H=4, C=6
MK-WPC64 ● Ø5.5, H=4, C=6

MK-NPK64 ● Ø4, H=4, C=6
MK-CPK64 ● Ø4.8, H=4, C=6
MK-WPK64 ● Ø5.5, H=4, C=6

MK-CPK81 ● Ø4.8, H=1, C=8
MK-WPK81 ● Ø5.5, H=1, C=8

MK-CPK82 ● Ø4.8, H=2, C=8
MK-WPK82 ● Ø5.5, H=2, C=8

MK-CPK83 ● Ø4.8, H=3, C=8
MK-WPK83 ● Ø5.5, H=3, C=8

MK-CPK84 ● Ø4.8, H=4, C=8
MK-WPK84 ● Ø5.5, H=4, C=8



Overdenture

LOCKiT Attachment System



Tightening
torque:
30Ncm

CONCAVE

NARROW

STANDARD

WIDE

H1

H2

H3

H4

H5



MN-LOC01



MN-LOC02




MN-LOC03



MN-LOC04



MN-LOC05




MD-LOC01




MD-LOC02




MD-LOC03




MD-LOC04



MD-LOC05



MW-LOC01



MW-LOC03



MW-LOC05



All LOCKiT products
shown in this catalog are
sold subject to individual
regulatory certification only.

LOCKiT System Components

Caps range



MM-BLA10
MM-RED10
MM-GRE10
MM-BLU10
MM-PIN10
MM-CLE10

Male processing package



MM-LOC02

Metal housing




MM-TIT10

Impression coping




MM-LOC1

Analog




MM-LOA40 - NP/SP
MM-LOA50 - WP

Abutment holder sleeve




MT-LOC11

Core tool



MT-LOC10

Ratchet driver



MT-RSD19 - Short
MT-RSD24 - Long

LOCKiT abutments are also sold in kits

MK-NLC1	● Ø3.75, H=1
MK-SLC1	● Ø3.75, H=1
MK-WLC1	● Ø4.5, H=1
MK-NLC2	● Ø3.75, H=2
MK-SLC2	● Ø3.75, H=2
MK-NLC3	● Ø3.75, H=3
MK-SLC3	● Ø3.75, H=3
MK-WLC3	● Ø4.5, H=3
MK-NLC4	● Ø3.75, H=4
MK-SLC4	● Ø3.75, H=4
MK-NLC5	● Ø3.75, H=5
MK-SLC5	● Ø3.75, H=5
MK-WLC5	● Ø4.5, H=5



Overdenture

OT-Equators

Tightening
torque:
30Ncm

NARROW

STANDARD


WIDE

H1




MK-NOE1

H2




MK-NOE2

H3




MK-NOE3

H4




MK-NOE4


H5




MK-NOE5




MK-SOE1




MK-SOE2



MK-SOE3




MK-SOE4




MK-SOE5



MK-WOE1



MK-WOE3



MK-WOE5


OT-Equator Components

Impression coping cap



OE-IC002

Analog



OE-RSM60

Metal housing and
transfer caps




OE-MH001
OE-RCB01
OE-RCY01
OE-RCP01
OE-RCW01
OE-RCV01

Handling tool




ET-IT002

Driver for 0.05 inch hex.



MT-RDS30 - Short
MT-RDL30 - Long

OT-Equators are sold in kits



Overdenture

Ball Attachment System

Tightening torque: 30Ncm

NARROW

STANDARD

WIDE


H1

H2

H3

H4


H5




MB-N2330




MB-N3330




MB-N4330




MB-N1375




MB-N2375




MB-N3375




MB-N4375



MB-N5375




MB-AN152




MB-AN252



MW-B1500




MW-B3500



MW-B5500

Ball Attachment System Components

Analog



MB-RS375

Titanium ball cap



MB-TTP10

Plastic mounting sleeve



MB-PR010

Gold ball cap



MB-GPS10

Plastic ball cap



MB-SF200

Cap spring




MB-TBS10

Cap spring



MB-TBS20

Standard plastic cap
for metal housing



MB-PPC10

Soft plastic cap for
metal housing




MB-PPR10

Plastic disc



MB-DB235

Double metal housing



MB-DMH10

Ball attachment
anchor key



MT-RB225

Prosthetic Screws and Drivers

i Tightening
torque:
20/30Ncm

Prosthetic screw



MN-S0160 - NP
MD-S0220 - SP

Medium prosthetic
screw



MD-S0222

Long prosthetic screw



MD-S0224

Mini prosthetic screw




MD-S0200

Gold prosthetic screw



MD-G0220

Short driver for 0.05
inch hex.



MT-RDS30

Long driver for 0.05
inch hex.



MT-RDL30

Extra-long driver for
0.05 inch hex.




MT-RDE30

Extra-short driver for
0.05 inch hex.




MT-RDM30

Short motor adapter
for 0.05 inch hex.



MT-SM005


Long motor adapter
for 0.05 inch hex.



MT-LM005


Prosthetic Screws and Keys for Angled Screw Channel

Ratchet key for angled
screw channel



MT-ELR10 - Long
MT-ESR10 - Short

Motor key for angled
screw channel




MT-ELM10 - Long
MT-ESM10 - Short

EZ-Base prosthetic
screw




MN-E0160 - NP
MD-E0220 - SP



A unique screw and driver tip
design enable extreme angulation
and a wide insertion path.

MK-0039
Prosthetic Tools Kit





All rights reserved. No part of this publication may be reproduced, transcribed, stored in an electronic retrieval system, translated into any language or computer language, or be transmitted in any form whatsoever, without the prior written consent of the publisher. Warning: MIS products referred to in this document should be used by licensed dentists only.

