

MIS | MGUIDE

Peace of Mind with Personalized Digital Dentistry

MIS[®]




MGUIDE is the industry's most advanced, yet simple solution for guided surgery. The 3D printed template is designed with an open-frame for maximum visibility, irrigation and accessibility from all angles without the need for removal. MGUIDE is a keyless system, designed for single handed procedures, eliminating the need for unnecessary tools. Raised flap surgery may also be more easily performed.

Increased visibility
through the open
frame design

Single-handed
procedure with the
keyless system

Access for irrigation
and anesthesia



Advantages



Accurate

The MGUIDE system incorporates a number of procedures designed to deliver accurate implant planning schematics and custom-designed guided surgery templates, for a more precise surgical procedure and predictable prosthetic results.



Less Chair Time

From start to finish, the MGUIDE system is shaped by precision state-of-the-art equipment, to help create an optimal plan and flawless template. This enables a smooth guided implant planning and placement procedure, cutting down on patient visits and ultimately saving valuable chair-time.



Smart Design

The advanced digitally designed template allows fully braced seating in the patient's mouth. The system also enables flexibility in performing the surgical procedure from an anatomical standpoint, including implant replacement as needed, without changing or redesigning the template.

MGUIDE Surgical Kits

The MGUIDE surgical kits are comprised of tools and instruments specially engineered to optimize the guided implant placement procedure.






The MISOFT program enables the creation of a highly accurate implant placement plan, taking all aspects of the patient anatomy into consideration, using a simple digital workflow. The MGUIDE template is then designed right in the software.

Flexible prosthetic-driven planning may be done either by the clinician, using our simplified MISOFT program, or through our worldwide network of MCENTERs, which provide highly professional technical support and online guidance.

MISOFT also acts as an online information hub connecting all software users; doctors, dental labs, periodontists, prosthodontists and the MCENTER in order to share cases and take part in demonstrations, discussions or consultations. The MISOFT virtual planning process is documented and stored on the Cloud for future reference.





MIS MCENTERS offer comprehensive, quality, digital dentistry services, including implant planning, MGUIDE template design and production, detailed surgical plans, as well as customized temporary solutions, all expertly engineered to deliver simplicity and convenience.



Find your local MCENTER
on the MIS Website

